



BEAR Micro Series Disposable Micro Vascular Clip

Features of BEAR's Disposable Vascular Clips

1. Single-use Reliability



BEAR's disposable clips are used only once – so clamping pressure is always exact and reliable.

2. Guaranteed Debris-Free, Corrosion-Proof, and Sterile

BEAR's clips never have to be scrubbed before use, and you will never find debris, rust, or corrosion.

3. Secure, Atraumatic Occlusion of Arteries and Veins

The jaws of each BEAR artery clip and vein clip have a unique hemispheric design that assures secure, non-slip occlusion.



Specification



The supervision of Prof. Susumu Tamai.

| Product code | Closing Power(gf) | Blood Vessel Diameter | Product size(mm) | | | |
|--------------|-------------------|-----------------------|------------------|-----|-----|-------|
| | | | A | B | C | D |
| TKF-1-15 | 15 | 0.2~1 | 7.0 | 4.0 | 2.3 | |
| TKF-1-30 | 30 | | 7.0 | 4.0 | 2.3 | |
| TKF-2-15 | 15 | | 7.0 | 4.0 | 2.3 | ~8.0 |
| TKF-2-30 | 30 | | 7.0 | 4.0 | 2.3 | ~8.0 |
| TKS-1-20 | 20 | MAX 1 | 10.0 | 4.5 | 2.5 | |
| TKS-1-40 | 40 | | 10.0 | 4.5 | 2.5 | |
| TKS-2-20 | 20 | | 10.0 | 4.5 | 2.5 | ~7.5 |
| TKS-2-40 | 40 | | 10.0 | 4.5 | 2.5 | ~7.5 |
| TKM-1-30 | 30 | 1~2 | 11.5 | 5.5 | 3.0 | |
| TKM-1-60 | 60 | | 11.5 | 5.5 | 3.0 | |
| TKM-2-30 | 30 | | 11.5 | 5.5 | 3.0 | ~11.5 |
| TKM-2-60 | 60 | | 11.5 | 5.5 | 3.0 | ~11.5 |
| TKL-1-60 | 60 | 2~4 | 15.0 | 6.5 | 3.5 | |
| TKL-1-120 | 120 | | 15.0 | 6.5 | 3.5 | |
| TKL-2-60 | 60 | | 15.0 | 6.5 | 3.5 | ~15.0 |
| TKL-2-120 | 120 | | 15.0 | 6.5 | 3.5 | ~15.0 |

How to select a product code

Body: polycarbonate

Spring: stainless steel

- ① Select a blood vessel diameter.
- ② Select a closing power, either for vein or artery.
- ③ Select a either single or double.



BEAR Medic Corporation

Tokyo Branch: Yushima Bear Bldg., 2-31-24,
Yushima, Bunkyo-ku, Tokyo, Japan 113-0034
Telephone: + 81 3 3818 4041 Fax: +81 3 3818 4042
URL: <http://www.bearmedic.co.jp/>