



BP Spatula
for plastic surgery

CC Reverse cutting
for plastic surgery

Needle type	Needle length(mm)	Suture type	USP
BP	6 ~ 15	nylon monofilament	8-0 ~ 3-0
	8 · 11 · 13	polypropylene	7-0 ~ 4-0
CC	8 · 11 · 13	nylon monofilament (black only)	7-0 ~ 5-0
		polypropylene	
E	6 ~ 24	nylon monofilament	8-0 ~ 2-0
		silk	7-0 ~ 4-0
		polypropylene	

E Reverse cutting
for plastic and orthopedic surgery

SP Spatula
for general surgery

(Reverse cutting for over 15mm length)

S Reverse cutting
for general surgery

Needle type	Needle length(mm)	Suture type	USP
SP	8 ~ 30	nylon monofilament	7-0 ~ 1-0
		silk	7-0 ~ 2号
		braided nylon	6-0 ~ 1-0
		tetoron	6-0 ~ 2-0
		polyester	6-0 ~ 2号
S	11 ~ 30	nylon monofilament	7-0 ~ 3-0



T Taper point
general taper point and microsurgery

ST S point
for cardiac surgery

Needle type	Needle length(mm)	Suture type	USP
T	3 ~ 40	nylon monofilament (black only for microsurgery)	11-0 ~ 1-0
		polyester	10-0 ~ 2号
		silk	6-0 ~ 3-0
ST	8 ~ 21	polypropylene	7-0 ~ 2-0

VT Tapercut
For cardiac surgery and microsurgery

T(M)WT Taper point
For double eyelid surgery


Needle type	Needle length(mm)	Suture type	USP
VT	4 ~ 35	nylon monofilament (black only)	11-0 ~ 8-0
		polypropylene	10-0 ~ 1-0
		polyester	7-0 ~ 1-0
T(M)WT	17 · 19 · 21	nylon monofilament (black and blue)	8-0 · 7-0



Tokyo Branch: 〒113-0034 Yushima Bear Bldg.,2-31-24,Yushima,Bunkyo-ku,Tokyo,Japan

Telephone:+81-3-3818-4041 Fax:+81-3-3818-4042

URL: <http://www.bearmedic.co.jp/>



F-20130225